JW Hospital Equipment

# SCRUB SILLIPURPOSE



MULTIPURPOSE

# AUTOMATIC SCRUB STATION

CHS-MATIC I, II, III









CHS-MATIC II CHS-MATIC I

CHS-MATIC Series are complete equipment with a special micro-filtering system, heating system for hot water and ultra-violet sterilizer. An ideal Scrub-Station which can be operated either manually or automatically.

AUTOMATIC FUNCTION OF CHS-MATIC SERIES ORIGINATES FROM THE FOLLOWING DATA CONTROL PANEL TO WHICH AN I.C ELECTRONIC CIRCUIT IS ATTACHED ENERGIZED AUTOMATICALLY ALL THE FUNCTIONINGS OF SCRUB-STATION.

### **Scope of Application**

- 1. Operating theater for general surgeries.
- 2. OB/GY theater/Delivery rooms(bathing infants)
- 3. Operating theater for orthopedic surgeries.
- 4. Operating theater for neuro-surgeries.
- 5. Operating theater for ENT surgeries.
- 6. Other places in need of sterilized water.

### **Features**

- 1. All automatic deluxe scrub station.
- 2. Sterilized water (99.6%) supply two U.V. water sterilizing lamps.
- 3. Hot water supply for four seasons.
- 4. Thorough filtering prefilter and main filter.
- 5. Easy operation by knee water and soap.
- 6. Water pressure control device.
- 7. Safe-dual temperature sensors for safety.
- 8. Rust-proof materials-all stainless steel (SUS304) and copper tubes.
- 9. Charming design.

### **Basic Operations and Function**

- 1. Heating
  - a. Hot water supply in winter.
  - b. Fast warming by one heater. (2KW)
  - c. Two safety thermo-switches to controll the in-creasing pressure of the warming tank.
  - d. Constant water temperature by safety thermo-switches.
- 2. Sterilization
  - a. The water can be sterilized to 99.6% (0.6T/H) by U.V. water sterilizing lamps.
  - b. The U.V. lamps are made of dual quart tubes for electrical safety.
  - c. Lamp use time can be read as life-timer is attached.
- 3. Filtering
  - a. Filtering to  $10\,\mu\text{m}$  by pre-filter.
  - b. Filtering to  $5\mu m$  by main-filter.
  - c. Thorough filtering by dual-filters.
  - d. Easy filter replacement.

## RESULT OF TEST PERFORMANCE SHOWS THAT THE STERILIZATION CAPABILITY IS AN AVERAGE 99.6%.

Scrub-Station's sterilization performance test by the use of ultraviolet sterilization method.

### 1. method


Germs Used for Tests	E.Coil, P.Aeruginasa, B.Subtilis.
Test Method	Numbers of germs are measured after diluted them into the refined water, while each germ of testing are collected more than 3 times from BHI Agarmodium under the state of 37°C and 2.4 Hrs.  Samples of each germ passed through S.S. unit in the amount of 5 liter/mm and 10 liter/mm are collected into the sterilized vessel of each and conducted testing in acordance with the rule of membrance filtering method as prescribed in the Korean Medicine General Testing Procedures.


### 2. Results

Germs, tested	Numbers of Germs before Passing thru S.S.	Passed thru	Percentage of sterilized(%)	Residuals
E.Coli	3.75×10⁵/mℓ	<b>5/</b> mℓ	99.99	none
P.Aeruginasa	3.75×10 <sup>5</sup> /mℓ	0.4/mℓ	100	none
B.Subtilis	1.25 × 10 ⁴/mℓ	1.5 × 10/mℓ	98.8	none

- 3. Findings: Average Sterilization Rate: 99.6%
  - · Superior in sterilization performance ability

### **DIMENSIONS**






No.	Part' s Name	Size	Q' ty	Remarks
1	Supply Pipe	15A	1	Hot Water
2	Supply Pipe	15A	1	Cold Water
3	Drain Pipe	65A	1	
4	Consent	110/220V	1	2.5Kw

### **MATTERS TO BE NOTED**

- 1. Please attach valve to supply pipe.
- 2. Capacity of breaker switch which is applied to the main switch is 30A
- 3. According to devices, for 3 persons use, 2,000mm 2 persons use, 1,400mm
  - 1 persons use, 1,000mm
- 4. Power source outlet of both for 3 persons use, 2 person's and 1 person's use are same.

### **FILTRATION**



# CONTROL PANEL FOR AUTOMATIC FUNCTIONS

### **PART'S NAME**

- ① U.V.R. Lamp operation indicator
- ② Water timing lamp
- ③ Water time switch
- (4) Heater Lamp
- ⑤ Heater on/off switch
- ⑤ Temperature control switch
- ① U.V.R. Lamp operation indicator
- Power lamp
- Power switch



# BRUSH DISPENSER CHS-BRS-D

### **SPECIFICATIONS**

Overall Size: W50mm×L120mm×H485mm Capacity: Brush 10ea, wall type

Material : 1.0t stainless steel plate

### **OPTIONAL ACCESSORIES**

Hand Scrub Brush



### **SPECIFICATION**

Division	Model	CHS-MATIC I	CHS-MATIC II	CHS-MATICⅢ	
$Size(W \times L \times H)$	mm	600 × 1,000 × 1,160	600×1,400×1,160	600 × 2,000 × 1,160	
Soap Pump	ea	1(12V-18W)	1(12V-18W)	3(12V-18W)	
Treatment		Germs removal by filtration			
Heater		2Kw×1ea			
Water	°C	0~60			
Filter	ea	Pre-Fiter 1, Main Fiter 1			
Sterilizer		Ultra Violet Ray Flow-Sterilizing System(15W × 2ea)			
Water Supply	kg/cm²	Min. 1.0-Max.3.0			
Power Source	V	AC220V, 50/60Hz			
Power Consumption	Kw	2.5			
Tank Capacity	l	Approx. 30			

OPTIONAL ACCESSORIES

Mirror: 1 set

The specifications are subject to change without prior notice for improvement



### Manufactured by

JW Tower, 2477 Nambusunhwan-ro, Seocho-gu, Seoul, Korea

**TEL** +82-2-840-6649 **FAX** +82-2-852-0565

**E-mail:** export@jw-holdings.co.kr www.jw-bioscience.co.kr